

PROTOCOL

1. Before beginning the treatment, make sure you have your [SAIAN Spray](#) and [SAIAN Cryo-Globes](#) in the refrigerator/freezer. The finishing cooling touch is very important at the end of this treatment.

2. Cleanse the face with [SAIAN Replenishing Cleanser](#). This milky cleanser's gentle and natural formula removes surface impurities, leaving skin moisturized and supple. Aloe Vera and chamomile are hydrating, softening, calming, as well as anti-microbial and anti-inflammatory. This product is perfect for sensitive skin. (You can also choose [SAIAN Strawberry Cleanser](#) or [SAIAN Botanical Lipid Cleanser](#))

3. Apply numbing cream to treatment area and leave on for 30 minutes. Remove numbing cream thoroughly.

4. Apply antiseptic wipe to treatment area.

5. Apply [1/2 SAIAN Hyaluronic Acid Ampoule](#) to the entire treatment area. The rest of the ampoule will be used later in the treatment.

6. Begin treatment by dividing the face into five parts: forehead, left and right cheeks, nose, upper lip and chin.

7. Perform microneedling according to your device's manufacturer's instructions. Work on all 5 sections in the same manner until the entire face has been treated. After treating each section, wipe any bleeding.

8. After the entire face has been treated, take the [SAIAN Pure Hyaluronic Acid Spray](#) out of the refrigerator, and spray the face, neck and décolleté. T

9. Cover the face with [SAIAN Restoring Mask](#)
or

Use our freeze dried [SAIAN Pure Collagen sheet Mask](#). Fold the mask in half over the client's face, mark and cut the mask around the nose area. Apply the dry collagen mask immediately over the face and spray with [SAIAN Pure Collagen Elastin Spray](#), molding gently to the face and neck with sponges. Leave mask on the face for no less than 20 minutes under steam, **NO OZON!** Remove the mask

10. Apply your choice of [SAIAN Creme](#) to the entire face and [SAIAN Clinical Strength Eye Serum](#) around the eyes.

11. Spray with [SAIAN Sunscreen Mist](#). This is an oil-free light sunscreen spray. Chemical sun blocks are never recommended after invasive treatments, and this light Titanium dioxide based sunscreen is the perfect alternative.

12. Roll over the entire area with [SAIAN Cryo Globes](#).

13: Recommend the client to use [SAIAN Natural Clinical Skin Care](#) at home to keep skin healthy and firm.

