

SAIAN®

NATURAL CLINICAL

ULTRASONIC
CAVITATION
CELLULITE
TREATMENT



ULTRASONIC CAVITATION CELLULITE TREATMENT 45 MIN

Using leading technology and our newest machine, SAIAN Mini Cavi Pro - ultrasonic cavitation, fat cells convert into liquid which then naturally can drain through the body's natural filtration system. This machine is ideal for use on the thighs, abdomen, love handles and buttocks. Generally, this treatment takes about 30 minutes to one hour, and 72 hours must pass between each treatment so that the body can eliminate the fat. We suggest at least 10-12 treatments for best results.

This ergonomic wireless handset is used in combination with Saian Bio-Contour Cellulite Slimming Gel, and will become your favorite tool!

The handset emits a low level of ultrasonic waves composed of the compression and expansion impulses traveling in high speed cycles. Ultrasound frequencies causing thousands of microscopic air bubbles to surround the lipocytes/fat cell membranes and through vibrations, creating bubbles with enough pressure to cause the breakdown of fat cells immediately. The content of the digested lipocytes then move through the lymphatic system to the liver and metabolizes in a safe manner and is excreted from the body, assisted by water intake both on the day of treatment and on subsequent days. The treatment works best at removing superficial fat deposits and fat deposits that are non-volumetric - you will absolutely love the result!

SAIAN Bio Contour Slimming Cellulite Gel is water-based and fragrance free, and is perfect to use with radio frequency equipment, cavitation machine, microcurrent, galvanic, and ultrasonic. It is also ideal for body wraps, thigh and abdomen wraps and arm wraps. This gel's revitalizing properties make it an obvious choice when it comes to natural slimming therapy.

PRODUCTS

- GLOVES
- Drybrushes
- Bio Contour Cellulite Slimming Gel
- Mini Cavi Pro machine

APPLICATION

1. Perform dry brushing on the area where the treatment will be done
2. Use the SAIAN **Mini Cavi Pro** machine with RED photon light and **ultrasonic cavitation**, level HIGH
3. Apply SAIAN **Bio Contour Cellulite Slimming Gel** to the machine head and to the area you will be working on.
4. Massage the area with circular motions towards the lymphatic drainage flow. 5 minutes per area. Massage remaining product into the skin with your hands.
5. Inform the patient to drink plenty of water in the coming days. 1 liter of water / 30 kg.

