

Our small, portable machine is easy to use, and as strong as large, stationary equipment

Our hand-held portable high frequency facial machine comes with 4 glass electrodes. When applied to the surface of the skin or scalp, a mild high frequency current passes through the neon or argon gas filled glass electrode causing it to light up with a calming orange-red or warm purple glow. The result is a pleasant, therapeutic tingling sensation. High frequency current causes an infusion of enriched oxygen molecules into the skin thereby cleansing, massaging, stimulating blood flow, re-texturizing and aiding in lymphatic drainage. The result is clearer skin.

The application of high frequency current promotes a cleansing antibacterial action that helps kill germs and bacteria which cause acne. The skin is left feeling instantly energized, rosy, firm, refreshed and noticeably softer – even after just one treatment.

ACNE: Treat blemished area(s) by lightly whisking the glass bulb tip of your high frequency wand across the surface of the skin using circular motions (treat enlarged pore areas in the same manner).

Concentrate on larger blemishes and the areas most affected for at least 30-40 seconds each using a gentle dabbing motion. Total treatment time is 3-5 minutes.

THINNING HAIR: Make sure scalp area and hair are clean. Set the intensity dial on your high frequency device to the highest setting. Touch the glass bulb tip of your high frequency device to the scalp and treat affected areas in a light circular or zig-zag motion for the recommended treatment time. For maximum effectiveness, apply hair re-growth or thickening product and gently massage into the scalp for 1-2 minutes. ALWAYS apply scalp treatment products FOLLOWING high frequency application.

CONTRAINDICATIONS:

Pregnancy; hypersensitive skin; diabetes; epilepsy; enlarged gland and thyroid problem; infectious, contagious, or acute condition/disease; with/had pacemaker or cardiovascular and/or cerebrovascular disorders; with metallic plate or anything metallic embedded anywhere in the body; with superficial cut, abrasion, injury, or scarred tissue; infectious or contagious skin conditions.

FAQ - How to use this machine correctly?

1. In order to avoid electrical shock during treatment, clients must not wear any jewelry or come in contact with any metal objects.

2. Do not use alcohol-based lotion together with hi-frequency treatment. If it is necessary to use an alcohol-based lotion, first treat the desired area with hi-frequency function function before applying the lotion.

3. Apply test the hi-frequency intensity before performing the procedure on the client.

4. Never use the hi-frequency function on the eyes, mouth, or any area with major and important gland. It is very important and the sole responsibility of a licensed and experienced esthetician to know the proper and forbidden areas to perform hi-frequency.

5. When performing the high frequency procedures, DO NOT keep the high frequency electrode on the same skin surface area too long for it may burn the skin and cause scars and dark pigmentation.

